



Mach7 Enterprise Imaging Platform v11.9.1

IHE Integration Statement



Manufacturer:
Mach7 Technologies, Inc.
120 Kimball Avenue, Suite 210
South Burlington, VT 05403
USA
+1 888-87-MACH7 (Toll Free USA
only)
+1 802-861-7745 (Phone)
+1 802-861-7779 (Fax)



**European Authorized
Representative:**
Emergo Europe
Prinsessegracht 20
2514 AP
The Hague, Netherlands
+31.70.345.8570 (Phone)
+31.70.346.7299 (Fax)

Australian Sponsor
Emergo Australia
Level 20, Tower II
Darling Park
201 Sussex Street
Sydney NSW 2000
Australia

All rights reserved. No parts of this work may be reproduced in any form or by any means -graphic, electronic, or mechanical, including photocopying, recording, taping, or information storage and retrieval systems without the written permission of the publisher. Products that are referred to in this document may be either trademarks and/or registered trademarks of the respective owners. The publisher and the author make no claim to these trademarks. While every precaution has been taken in the preparation of this document, the publisher and the author assume no responsibility for errors or omissions, or for damages resulting from the use of information contained in this document or from the use of programs and source code that may accompany it. In no event will the publisher and the author be liable for any loss of profit or any other commercial damage caused or alleged to have been caused directly or indirectly by this document.



Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or medical institution.

Technical Support

For technical queries, write to support@mach7t.com or visit our support site:
<http://support.mach7t.com>

General Enquiries



See User Guides / Instructions for Use

For general enquiries, write to info@mach7t.com.

Contact MACH7 Technologies

Americas Sales & Service

Mailing Address:

PO Box 586

Burlington, VT 05402

Office Location:

120 Kimball Avenue, Suite 210

South Burlington, VT 05403

T: +1-888-87-MACH7 (Toll Free USA only)

T: +1-802-861-7745

F: +1-802-861-7779

Release History

Release Version	Release Date	Sections Affected	Description
11.7.2	September 2016	General Release	This is the general release of this statement for the Mach7 Enterprise Imaging Platform v11.7.2
11.8.3	November 2018	General Release	This is the general release of this statement for the Mach7 Enterprise Imaging Platform v11.8.3
11.8.4	August 2019	General Release	This is the general release of this statement for Mach7 Enterprise Imaging Platform v11.8.4
11.9.0	July 2020	Conformance Statements Links	Updated link to other conformance statements.
v11.9.1	January 2021	General Release	This is the general release of this statement for Mach7 Enterprise Imaging Platform v11.9.1

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Audit Trail and Node Authentication (ATNA)	Secure Node	none
Access to Radiology Information (ARI)	Image Manager	none
	Image Archive	none
	Report Repository	none
	Image Manager	none

Consistent Presentation of Images (CPI)	Image Archive	none
Consistent Time (CT)	Time Client	none
Cross-Enterprise Document Sharing for Imaging (XDS-I.b)	Document Source	none
	Document Consumer	none
	Document Repository	none
Cross-Enterprise Document Sharing for Imaging (XDS)	Document Repository	none
	Document Registry	none
Evidence Documents (ED)	Image Manager	none
	Image Archive	none
Image Object Change Management (IOCM)	Image Manager	none
	Image Archive	none
	Image Display	none
	Change Requester	none
Key Image Note (KIN)	Image Manager	none
	Image Archive	none
Mammography Image (MAMMO)	Image Manager	none
	Image Archive	none
Multiple Image Manager Archive (MIMA)*	Image Manager	none
	Image Archive	none
Nuclear Medicine Image (NMI)	Image Manager	none
	Image Archive	none

Patient Identifier Cross-Reference HL7 V3 (PIXV3)	Identity Consumer	none
Patient Information Reconciliation (PIR)*	Image Manager	none
	Image Archive	none
Portable Data for Imaging (PDI)	Image Display	none
	Portable Media Importer	none
	Portable Media Creator	none
Radiation Exposure Monitoring (REM)	Image Manager	none
	Image Archive	none
Scheduled Workflow (SWF)*	Image Manager	none
	Image Archive	none
Teaching File and Clinical Trials (TCE)	Export Manager	none
Link to Standards Conformance Statements for the Implementation		
HL7	https://www.mach7t.com/resources/conformance-statements/	
DICOM		

*Supported profiles which may require customizations for Modality Procedure Step transactions.