

# eUnity

---

by  MACH7

eUnity 7.6

Enterprise Viewing, Integration, and Collaboration Platform

## IHE Integration Statement

Document ID: EUNITY-024

Revision number: 1.4

Revision date: February 5, 2025

# Table of contents

---

IHE integration statement .....	3
---------------------------------	---

# IHE integration statement

Product: eUnity

Vendor: Mach7 Technologies Canada Inc.

Integration profiles implemented	Actors implemented	Options implemented
Scheduled Workflow (SWF)	Image display	None
Consistent Presentation of Images (CPI)	Image display	None
Simple Image and Numeric Report (SINR)	Report reader	None
Access to Radiology Information (ARI)	Report reader	Multiple sources
	Image display	Multiple sources
Audit Trail and Node Authentication (ATNA)	Secure application	Encryption option
	Audit record repository	None

## Links to general information on IHE

North America	<a href="http://www.ihe.net">www.ihe.net</a>
Europe	<a href="http://www.ihe-europe.net">www.ihe-europe.net</a>
Japan	<a href="http://www.ihe-j.org/">http://www.ihe-j.org/</a>

## eUnity DICOM conformance

For the eUnity DICOM conformance statement, see <https://mach7t.com/collateral>.